

CROWN UNDER PARTIAL PROCEDURE

1. Take **first** impression *without the partial* in place. This is to be used for the removable die and accurate margins (can be a quadrant).
2. Take **second** impression (full mouth) *with the partial in place* (using the same type crown and bridge impression material used in the first impression).
3. Partial **must** remain in impression until poured in stone. Send the partial and impressions to lab for processing.

If the patient can not go without the partial for a few days, then block out the undercuts at the acrylic saddle areas with wax and pour the impression with a die stone. After bench set, carefully remove the partial from the model.

After the initial model and die work is completed in the laboratory, we will need the partial twice. Usually, patients can wait at the lab while we utilize the partial, generally the procedures take 30 minutes.

If you have any questions, please call.